

Available Wheelchair Lifts for Full-Size Ford Vans

The **Braun Vangater Series Wheelchair Lifts** are fully electric in operation, for both the fold/unfold and up/down cycles. Lift operation is controlled by the standard hand-held control, the on-lift controls, or the optional remote control. The Vangater Series is designed for installation in the side doors of domestic full-size vans. A model is also available for installation in minivans. With the compact tri-fold platform, the lift occupies minimal space inside the vehicle.



The **Braun UVL Series Wheelchair Lifts** are electric in operation for both the in/out cycle and the outboard roll stop. The up/down cycle is hydraulic in operation. All functions are controlled by the standard hand-held control or the optional remote control. The UVL Series is designed for installation at the side door of domestic full-size vans. Installation of the UVL requires modification to the van's drive shaft.



The **Braun Millennium Series Wheelchair Lifts** are the true workhorses of the BraunAbility line. Designed for installation in the side or rear doors of domestic full-size vans, three standard models with platform lengths up to 51" ensure there is a lift to fit your needs. With dual hydraulic lifting arms, you'll come to rely on the lift's strength and dependability day after day, year after year. The Millennium also has a side-entry cut option on the platform for easier loading/unloading in tight places (not available on the Century).



The **Braun Century Series Wheelchair Lift** is an economical solution to personal mobility. Like the Millennium Series, the Century Series lifts are fully hydraulic in operation, for both the fold/unfold and up/down cycles. Lift operation is controlled by the standard hand-held control, the on-lift controls, or the optional remote control.



The **VMI/Ricon Clearway Wheelchair Lift** offers an exclusive split-platform that folds vertically, allowing unobstructed access upon exit and entry whether or not the lift is deployed. The driver's line of sight is also unobstructed for optimum safety. Dual power cylinders ensure smooth, reliable operation, even for heavy wheelchairs and scooters. The strong all-steel frame and powerful hydraulic pump allow quick, reliable operation of the lift.



The **VMI/Ricon KlearVue Platform Lift** offers a unique horizontal fold to allow an unobstructed view from inside or outside the vehicle without requiring a raised roof or lowered floor. No vehicle modification is needed thanks to advanced engineering that allows a longer 51" platform to fit in a shorter 48" door.



The **VMI/Ricon Reliant Wheelchair Lift** features a wide range of platform sizes and can accommodate a mobility device up to 51". With a welded all steel frame, a powerful hydraulic pump, and a simple electrical system for trouble-free operation, Reliant lifts are a solid, strong and proven mobility solution. It features patented technology that securely locks the platform in the stowed position, automatic inboard/outboard roll stops, a non-skid platform, and convenient handrails with a control switch for lowering and raising the lift.



The **VMI/Ricon Unilite Wheelchair Lift** features an internal swing door operator and sliding door operator, a powered rollstop with mechanical latch, and handrails with control switch. The Unilite also features a built-in manual back-up system and may be equipped with an optional remote control.



The **VMI/Ricon Slide-Away Platform Lift** features strength and space in one easy-to-use package. When in use, the sliding second tower travels toward the front of the vehicle to provide the strength and stability only available in a dual post lift. And when not in use, the sliding second tower travels towards the rear of the vehicle, leaving the side door unobstructed or able-bodied passengers and providing full use of the front passenger seat.



The **VMI Fiorella F500 Platform Lift** is the newest platform lift in the VMI line of mobility transportation solutions. Fiorella has manufactured wheelchair lifts for more than 15 years and distributes in 47 different countries. This Fiorella platform lift features a unique design that is attuned to the specific needs of disabled individuals. Customer experience is the primary focus behind the engineering and creativity that goes into each and every F500 platform mobility lift.

